

1) Product Overview and Details

- The iCare HOME2 is a rebound tonometer used for measuring intra-ocular pressure (IOP) in glaucoma patients.
- IOP measures have a diurnal rhythm and 75% of IOP elevations occur outside of normal clinic hours (Liu et al, 2020).
- Glaucoma often goes undiagnosed in 50% of patients due to the IOP variability, thus home IOP monitoring has become the new standard to (Ertel et al, 2021):
- 1. Diagnose and Prevent disease progression.
- Monitor medication adherence and treatment regimes. Monitor post surgical outcomes.

Technical Details

Setup





Each tonometer has a carrying case.

The iCare HOME2 is the newest version of the iCare product line of home tonometer's.

- iCare HOME2 tonometer
 - Probe applicator USB cables for mobile and PC connection
 - USB stick with instructions for use
 - 4 x AA 1.5 V alkaline batteries






iCare, 2023


Module 3 - Portfolio Assignment: Digital/Virtual Health Product Infographic

Specialized Graduate Certificate in Health Informatics Concentration: Digital Health

Adam M. Lang University of Denver University College October 1, 2023


Faculty: Matt Vestal, M.H.A. Director: Jeffrey Weide, DBA Dean: Michael J. McGuire, MLS






• Unlike other Tonometer's, you can: sit, stand, or lye down.

User inserts sterile probe through the probe applicator.



to your face.




Patient Prompting



Measuring the IOP

Green ring notifies patient of optimal eye position.






adjust device to obtain optimal IOP measurements.



- iCare HOME2 positioned 90 degrees to face. User will press side buttons 6 times or hold.
- The green circle lights up with a "beep" for each of the 6 IOP measurements performed.
- Mobile Application: iCare PATIENT2




iCare PATIENT2



iPhone and Android app stores.

App available for download via




destination such as computer or cloud.

icare CLINIC # Pati



· User interface is friendly.



2) Data Collection and Processing Details

1. iCare PATIENT2 mobile application. **2. iCare CLINIC** cloud database and clinical application. 3. iCare EXPORT desktop application (not pictured below).

Data collected by the iCare HOME2 is transmitted to 3 platforms:

Patient and Provider can review the following data: IOP measurements for each eye in mmHg.

mobile app.


- Time of day measurements were performed. • Date measurements were performed.
- Quality of measurements based on device position to each eye. Two images **below** show data displayed in the **iCare PATIENT2**


1. Line graph for each eye, Right vs. Left. 2. Line graph with day/time vs. IOP measure. 3. Individual data points in a chart for each day, time, and each

eye.

• We can see the diurnal IOP measures displayed:

Aug 25 https://doi.org/10.1097/IJG.0000000000002298





8:59 PM • 13 Right




Chart type Diurnal Example of a "Diurnal" chart or graph within the iCare CLINIC cloud application. There are multiple graph and chart types that can be used to review IOP measures.

3) Product Accuracy and Effectiveness - Evidence Based Research


Time/Day Device used Quality of IOP measures

Dashboard to the <u>left</u> is the provider's iCare CLINIC

Providers can filter data points they want to

review and visualize based on:

cloud app.



Example of a line graph within the iCare CLINIC cloud app. Patient and Provider can review IOP measures for each eye, date and time.

Kratz et al. evaluated 135 eyes in 70 patients confirming the new iCare HOME2 tonometer is just as effective at measuring IOP as the gold standard "Goldmann applanation tonometer" or "GAT". • The graph below shows just how similar the IOP measures were between the 2 devices head-to-head: Mean IOP for the GAT was: 16.3 +/- 6.5 mmHg (range 3-56) • Mean IOP for the iCare HOME2 was: 16.5 +/- 7.3 mmHg (range 3-55)


Figure 1. Correlation of intraocular pressure (IOP, mmHg) measurements by Haag Streit AT900 GAT (IOPG) and iCare HOME2 tonometer (IOPi). Strong, significant correlation was found

Figure 1: Assaf Kratz, Rabia Zbidat, Rina Kishner, Michal Cohen, Walid Shalata, Ivan Goldberg, 2023. "Assessment of the iCare HOME2, a New IOP Self-Measurement Tonometer." J Glaucoma 2023

- Elyse J McGlumphy, MD, Aleksandra Mihailovic, ScM, Pradeep Y Ramulu, MD, PhD, Thomas V Johnson, MD, PhD, between paired IOP measurements from the two instruments. (IOP- intra ocular pressure, GAT-2021. "Home self-tonometry trials compared to clinic tonometry in patients with glaucoma." Ophthalmol Glaucoma 4, no. 6 (Nov-Dec 2021): 569-580, https://doi.org/10.1016/j.ogla.2021.03.017 Goldmann applanation tonometer)
 - McGlumphy et al. examined 107 eyes in 61 glaucoma patients. The maximum IOP occurred outside of the 8am to 5pm clinic hours on ~50% of the days and between 4:30am to 8am on 24% of days. The maximum daily mean IOP significantly exceeded those seen in the clinic in ~50% of patients. The graph below shows a patient that had recurring IOP elevations averaging 25-40 mmHg every morning between 6am-7am yet their IOP measures during clinic hours were in the low teens.

This study validated:


McGlumphy et al.



Efficacy of using the iCare HOME to monitor IOP spikes that commonly occur outside of clinic.

to years at a time if they never have an IOP spike during clinic hours.

The accuracy of using the iCare HOME to diagnose patients with glaucoma who will go undiagnosed for months



https://doi.org/10.2147/OPTH.S284844 Liu et al. published case reports of using the iCare HOME tonometer for pre and post operative IOP monitoring in glaucoma patients. The pre and post operative IOP graphs for "Case 4" to the left, shows a 61 yr. old male with consistent IOP elevations outside of clinic hours at midnight for both eyes, diagnostic of open angle glaucoma.

The authors concluded from their cohort that the iCare HOME tonometer shows accuracy and efficacy in the pre and post operative setting for monitoring IOP measures in glaucoma patients.

Patient had surgical intervention only in the left eye (OS, red line) and the iCare HOME tonometer demonstrated the improved outcome from the procedure with the red line IOP measures returning to standard baseline during diurnal

1. <u>iCare CLINIC</u> – cloud browser-based program and database for storing and viewing all IOP measurements, transferring data to external systems such as EHRs and for clinical decision making.

There are 3 platforms involved in Data transmission:

4) Flow of Data Transmission

iCare EXPORT – desktop application used to transport IOP data to the iCare CLINIC database or to an external system such as an EHR or computer. It can also display IOP measurement data. iCare PATIENT2 - tablet/mobile application that receives IOP data from the iCare HOME2 tonometer, transfers IOP data to the

iCare CLINIC database or to an external system of choice.

monitoring while the right eye continued to show IOP spikes.

- Diagram of Data Flow Transmission (iCare, 2023) • **Top left:** data transferred via Bluetooth or USB to mobile application
- Micro USB USB-A the iCare PATIENT2 mobile app. • Bottom: Direct transfer of data from the iCare HOME2

USB-C

Diagram of USB data transmission (iCare, 2023) • Top: Direct transfer of data from the iCare HOME2 to

- Page 18 → Day 1
- Review of Characteristics and Clinical Utility." Clin Ophthalmol 2020, no. 14 (Nov 23): 4031-4045,
- iCare PATIENT2. <u>Top right:</u> data transferred via Bluetooth or USB to iCare EXPORT desktop PC application.
- to the iCare EXPORT. **Bottom:** Data is synced and stored in the **Care CLINIC** database application where it can also be reviewed in a browser.